

MATTERS ARISING

Open Access



Is qPCR always the most sensitive method for malaria diagnostic quality surveillance?

Cristian Koepfli^{1*}

Abstract

In many studies to evaluate the quality of malaria diagnosis, microscopy or rapid diagnostic tests (RDT) are compared to PCR. Depending on the method for sample collection and storage (whole blood or dried blood spot), volume of blood used for extraction, volume of DNA used as PCR template, and choice of PCR target (single vs. multi-copy gene), the limit of detection (LOD) of PCR might not exceed the LOD of expert microscopy or RDT. One should not assume that PCR always detects the highest number of infections.

Background and discussion

In a recent study, Abebe and colleagues found a large number of patients diagnosed to carry *Plasmodium falciparum*/*Plasmodium vivax* mixed infections by microscopy or rapid diagnostic test (RDT) to be mono-infections by qPCR [1]. Among febrile patients in Ethiopia, only 20/34 mixed infection detected by RDT, and 18/68 mixed infection detected by expert microscopy were confirmed by qPCR [1]. The authors propose that RDT and expert microscopy yielded incorrect results, and called for a strengthening of routine malaria diagnostic methods. One needs, however, to question whether qPCR could have missed mixed infections.

The limit of detection (LOD) of molecular diagnostics is impacted by the type of input sample used, e.g., dried blood spot (DBS) or whole blood, blood volume analysed, DNA extraction protocol, and qPCR assay. The authors have provided the details needed to calculate a rough estimate of the LOD of their assay. They extracted DNA from one 6 mm-punch of DBS using the QIAamp DNA extraction kit and eluted the DNA into a

volume of 100 μ L. One 6 mm punch of DBS is expected to harbour the equivalent of approximately 15 μ L blood [2], thus the DNA was diluted approximately sevenfold during extraction (from 15 to 100 μ L). Holzschuh and Koepfli previously showed that only approximately 10% of DNA is recovered from a DBS during the lysis step, and that the QIAamp kit recovers approximately 25% of DNA [3]. Considering dilution and loss of DNA during extraction, the number of genomes per μ L extracted DNA was reduced by a factor of 280 compared to initial parasite density (sevenfold dilution * tenfold reduction of DNA during lysis of DBS * fourfold reduction during DNA clean up). The study authors used 5 μ L of extracted DNA as template for qPCR and amplified a mitochondrial gene [4], which is expected to be present in around 10 copies per parasite [5]. Let's assume an infection with an initial density of 10 parasites/ μ L, resulting in 100 mitochondrial genomes per μ L blood. After extraction, 0.35 (100/280) mitochondrial genomes are expected per μ L DNA. The 5 μ L DNA used as template correspond to less than 2 mitochondrial genomes (5 * 0.35) added to the qPCR. It is difficult to determine the minimal number of templates that a PCR is able to amplify, and thus whether the sample in this example would reliably result in a positive PCR result. Many replicates with well-quantified template densities around the limit of amplification would need to be run to determine that. While some PCRs are able to yield

*Correspondence:

Cristian Koepfli
ckoepfli@nd.edu

¹ Department of Biological Sciences and Eck Institute for Global Health, University of Notre Dame, Notre Dame, USA



a positive result from as little as one template, others require higher template densities.

Even with the actual efficacy of DNA extraction of the current study likely to differ from published estimates due to small differences in the protocols, it is possible that the LOD of the qPCR was well above 10 parasites/ μ L blood. Both expert microscopy [6, 7] and *P. falciparum* RDTs [8, 9] might detect parasites at densities as low as 10–50 parasites/ μ L. Expert microscopy and RDT might thus reach similar, or better LOD than qPCR (of note, Abebe and colleagues used the SD Bioline Malaria Ag Pf/Pv RDT, which is less sensitive than the latest generation of RDTs).

Challenges in diagnosis are even more pronounced among asymptomatic infections that are typically of low density [10]. Co-infections of *P. falciparum* and *P. vivax* might result in modulation of parasite densities [11, 12], thus adding to the discrepancies in the detection of co-infections by different diagnostic tools by Abebe and colleagues. When comparing microscopy and PCR to RDT data, further complexities arise as the concentration of antigens detected by RDTs do not always correlate with parasite densities [9]. Antigens can persist after parasite clearance [13], which could also explain some of the mixed infection detected by RDT where only one species was detected by PCR.

In summary, Abebe and colleagues highlighted the importance of considering mixed-species infections when the quality of malaria diagnosis is assessed. As for other studies that extracted DNA from a comparably small volume of blood collected as DBS [14], it is important to remember that well-trained expert microscopists might reach a LOD similar or better than PCR. A calculation of the number of target molecules that are actually present in the aliquot that is used for the amplification reaction should be provided whenever an amplification protocol is presented. This includes a calculation of the corresponding volume of blood used as template for the reaction, considering dilutions made during extraction. Especially when a discussion of the limit of detection is presented, efforts should be made to quantify the proportion of genomes that are lost during extraction, for example by comparing densities determined by expert microscopy to target gene quantification by qPCR. Further, the number of copies per genome of the targeted gene needs to be reported. In this way, it will be possible to assess accurately the real limit of detection of any particular method. qPCR should not automatically be considered the most sensitive method to detect infections.

Acknowledgements

Not applicable.

Author contributions

CK conceived and wrote this manuscript.

Funding

No funding was received for this work.

Availability of data and materials

Not applicable.

Declarations

Ethics approval and consent to participate

Not applicable.

Consent for publication

Not applicable.

Competing interests

The author declares that no competing interests exist.

Received: 25 August 2023 Accepted: 13 December 2023

Published online: 15 December 2023

References

1. Abebe A, Ménard D, Dugassa S, Assefa A, Juliano JJ, Lo E, et al. Significant number of *Plasmodium vivax* mono-infections by PCR misidentified as mixed Infections (*P. vivax/P. falciparum*) by microscopy and rapid diagnostic tests: malaria diagnostic challenges in Ethiopia. *Malar J.* 2023;22:201.
2. WHO. Preparation of blood spots on filter paper. Geneva: World Health Organization; 2016.
3. Holzschuh A, Koepfli C. Tenfold difference in DNA recovery rate: systematic comparison of whole blood vs. dried blood spot sample collection for malaria molecular surveillance. *Malar J.* 2022;21:88.
4. Canier L, Khim N, Kim S, Sluydts V, Heng S, Douring D, et al. An innovative tool for moving malaria PCR detection of parasite reservoir into the field. *Malar J.* 2013;12: 405.
5. Britton S, Cheng Q, Grigg MJ, Poole CB, Pasay C, William T, et al. Sensitive detection of *Plasmodium vivax* using a high-throughput, colourimetric loop mediated isothermal amplification (HtLAMP) platform: a potential novel tool for malaria elimination. *PLoS Negl Trop Dis.* 2016;10: e0004443.
6. WHO. Basic malaria microscopy. Geneva: World Health Organization; 2010.
7. Payne D. Use and limitations of light microscopy for diagnosing malaria at the primary health care level. *Bull World Health Organ.* 1988;66:621–6.
8. Niyukuri D, Sinzinkayo D, Troth EV, Oduma CO, Barengayabo M, Ndereyimana M, et al. Performance of highly sensitive and conventional rapid diagnostic tests for clinical and subclinical *Plasmodium falciparum* Infections, and hrp2/3 deletion status in Burundi. *PLoS Glob Public Health.* 2022;2: e0000828.
9. Das S, Jang IK, Barney B, Peck R, Rek JC, Arinaitwe E, et al. Performance of a high-sensitivity rapid diagnostic test for *Plasmodium falciparum* malaria in asymptomatic individuals from Uganda and Myanmar and naive human challenge infections. *Am J Trop Med Hyg.* 2017;97:1540–50.
10. Bousema T, Okell L, Felger I, Drakeley C. Asymptomatic malaria infections: detectability, transmissibility and public health relevance. *Nat Rev Microbiol.* 2014;12:833–40.
11. Holzschuh A, Gruenberg M, Hofmann NE, Wampfler R, Kiniboro B, Robinson LJ, et al. Co-infection of the four major *Plasmodium* species: effects on densities and gametocyte carriage. *PLoS Negl Trop Dis.* 2022;16: e0010760.
12. Koepfli C, Nguitragool W, de Almeida ACG, Kuehn A, Waltmann A, Kattenberg E, et al. Identification of the asymptomatic *Plasmodium falciparum* and *Plasmodium vivax* gametocyte reservoir under different transmission intensities. *PLoS Negl Trop Dis.* 2021;15: e0009672.
13. Jang IK, Tyler A, Lyman C, Kahn M, Kalnoky M, Rek JC, et al. Simultaneous quantification of *Plasmodium* antigens and host factor C-reactive protein

in asymptomatic individuals with confirmed malaria by use of a novel multiplex immunoassay. *J Clin Microbiol.* 2019;57:e00948-18.

14. Russell TL, Grignard L, Apairamo A, Kama N, Bobogare A, Drakeley C, et al. Getting to zero: micro-foci of malaria in the Solomon Islands requires stratified control. *Malar J.* 2021;20:248.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Ready to submit your research? Choose BMC and benefit from:

- fast, convenient online submission
- thorough peer review by experienced researchers in your field
- rapid publication on acceptance
- support for research data, including large and complex data types
- gold Open Access which fosters wider collaboration and increased citations
- maximum visibility for your research: over 100M website views per year

At BMC, research is always in progress.

Learn more biomedcentral.com/submissions

