

METHODOLOGY

Open Access



Streamlined, PCR-based testing for *pfhrp2*- and *pfhrp3*-negative *Plasmodium falciparum*

Jonathan B. Parr^{1*†}, Olivia Anderson^{2†}, Jonathan J. Juliano^{1,2,3} and Steven R. Meshnick^{1,2}

Abstract

Background: Rapid diagnostic tests (RDTs) that detect histidine-rich protein 2 (PfHRP2) are used throughout Africa for the diagnosis of *Plasmodium falciparum* malaria. However, recent reports indicate that parasites lacking the *pfhrp2* and/or histidine-rich protein 3 (*pfhrp3*) genes, which produce antigens detected by these RDTs, are common in select regions of South America, Asia, and Africa. Proving the absence of a gene is challenging, and multiple PCR assays targeting these genes have been described. A detailed characterization and comparison of published assays is needed to facilitate robust and streamlined testing approaches.

Results: Among six *pfhrp2* and *pfhrp3* PCR assays tested, the lower limit of detection ranged from 0.01 pg/μL to 0.1 ng/μL of *P. falciparum* 3D7 strain DNA, or approximately 0.4–4000 parasite genomes/μL. By lowering the elongation temperature to 60 °C, a tenfold improvement in the limit of detection and/or darker bands for all exon 1 targets and for the first-round reaction of a single exon 2 target was achieved. Additionally, assays targeting exon 1 of either gene yielded spurious amplification of the paralogous gene. Using these data, an optimized testing algorithm for the detection of *pfhrp2*- and *pfhrp3*-negative *P. falciparum* is proposed.

Conclusions: Surveillance of *pfhrp2*- and *pfhrp3*-negative *P. falciparum* requires careful laboratory workflows. PCR-based testing methods coupled with microscopy and/or antigen testing serve as useful tools to support policy development. Standardized approaches to the detection of *pfhrp2*- and *pfhrp3*-negative *P. falciparum* should inform efforts to define the impact of these parasites.

Keywords: Rapid diagnostic tests, False-negative, Diagnostic resistance, Histidine-rich protein, hrp2, hrp3, RDT, Deletion, Malaria, *Plasmodium falciparum*

Background

Diagnostic testing is a core component of recent malaria control efforts. In Africa, where the majority of deaths due to malaria occur, rapid diagnostic tests (RDTs) are the most commonly employed malaria diagnostic strategy, accounting for 74% of diagnostic testing among suspected malaria cases [1]. The most commonly used RDTs in Africa rely upon detection of PfHRP2, a *Plasmodium falciparum*-specific antigen expressed by the histidine-rich protein 2 (*pfhrp2*) gene. However, recent reports from select locations in South America, Asia, and Africa

of *P. falciparum* parasites lacking *pfhrp2* and/or the histidine-rich protein 3 (*pfhrp3*) gene, which produces an antigen that cross reacts with some PfHRP2-based RDTs, raise concerns about the effectiveness of PfHRP2-based RDTs in affected regions [2–18]. In response, the World Health Organization (WHO) has prioritized efforts to address parasites with deletions of the *pfhrp2* and/or *pfhrp3* (*pfhrp2/3*) genes [19, 20].

The methods required to identify and confirm *pfhrp2/3* gene deletions are challenging, due to the difficulty of proving the absence of a gene. While PCR assays that target *pfhrp2/3* are expected to yield negative results when applied to parasites lacking the gene(s), PCR failure can occur for other reasons. Testing of parasites with intact *pfhrp2/3* genes may yield false-negative results due to DNA concentrations below the assay's limit of detection,

*Correspondence: jonathan_parr@med.unc.edu

†Jonathan B. Parr and Olivia Anderson are co-first authors

¹ Division of Infectious Diseases, Department of Medicine, University of North Carolina, 130 Mason Farm Rd, Chapel Hill, NC 27599, USA

Full list of author information is available at the end of the article

poor quality DNA, variable reagent performance, or other factors.

Cheng et al. published useful guidelines to standardize the reporting of *pfhrp2/3* gene deletions [21]. However, the specific methods employed for the detection and confirmation of deletions continue to vary between laboratories, and recent evidence suggests that atypical elongation temperatures may improve amplification of AT-rich regions of both genes [14, 22]. This manuscript seeks to address these issues by comparing the performance of published PCR assays for *pfhrp2* and *pfhrp3*, exploring the impact of reduced elongation temperatures on assay sensitivity, assessing assay specificity, and describing a streamlined testing algorithm.

Methods

The performance characteristics of six published PCR assays, including four designed to amplify *pfhrp2* and two designed for *pfhrp3*, were compared [3, 5, 9, 23]. After determining the optimal annealing temperatures for each assay, we assessed their lower limits of detection (LOD) using DNA extracted from cultured *P. falciparum* 3D7 strain parasites. DNA was quantified using the Qubit 2.0 instrument with dsDNA high sensitivity reagents (ThermoFisher Scientific, Waltham, MA) and serially diluted in nuclease-free water to achieve concentrations ranging from 10^{-1} to 10^{-7} ng/ μ L (seven tenfold dilutions). For each dilution, the assay was performed in triplicate, using different elongation temperatures of 60, 65, and 72 °C. PCR assays were performed on Mastercycler thermocyclers (model AG 22331; Eppendorf, Hamburg, Germany) using 25 μ L reaction volumes containing 12.5 μ L HotStarTaq Master Mix (Qiagen, Venlo Netherlands), 200–400 nM primers synthesized by Eurofins Genomics (Louisville, KY) with salt-free purification, nuclease-free water, and 3 μ L of DNA template (Table 1). For nested reactions, first-round product was diluted 100-fold in nuclease-free water prior to second-round amplification. PCR products were visualized using electrophoresis with 1% agarose gels in TBE buffer (Tris/Borate/EDTA). Finally, LOD testing was repeated using optimized reaction conditions and serial dilutions of 3D7 DNA from a separate stock.

The specificity of the best performing assays, including those with targets spanning exon 1 and 2 (exon 1/2) and exon 2 alone, was then evaluated. Assays were performed using control DNA from *P. falciparum* Dd2 (MRA-150G) and HB3 (MRA-155G) strain parasites, which lack *pfhrp2* and *pfhrp3*, respectively. Control DNA was obtained from the Malaria Research and Reference Reagent Resource Center ([MR4], BEI Resources, Manassas, Virginia) and diluted to a concentration of 0.1 ng/ μ L after initial quantification using Qubit as above. For

assays that yielded an unexpected result using optimized reaction conditions (i.e. bands from a *pfhrp2* assay performed using *pfhrp2*-deleted Dd2 DNA or bands from a *pfhrp3* assays performed using *pfhrp3*-deleted HB3 DNA), amplicons were sequenced using Sanger sequencing at Eton Bioscience (Research Triangle Park, NC), and assays were repeated at the other elongation temperatures (60, 65, and/or 72 °C). For PCR products with multiple bands appreciated by gel electrophoresis, individual bands were excised, and DNA was extracted using the QIAquick Gel Extraction Kit (Qiagen, Hilden, Germany) before sequencing. Gel extraction was performed according to the manufacturer's instructions, with the exception of the final DNA elution step, in which we performed two separate elutions using 30 μ L aliquots of Buffer EB through the column, followed by a final elution step using the combined 60 μ L initial eluate to maximize DNA yield. Raw sequence reads were processed using Sequencher 5.4 (Gene Codes Corporation, Ann Arbor, MI), trimming bidirectional sequences based on confidence values and visual inspection of the chromatograms. We used EMBOSS Water for pairwise nucleotide alignments to the 3D7 (v3.0) *pfhrp2* and *pfhrp3* reference sequences. Sequence identification was based on sequence homology and alignment score, using default settings (DNAfull matrix, gap open penalty 10, gap extend penalty 0.5) [24].

Results

Assay performance

Reduced elongation temperatures improved the sensitivity of five of the six assays (Additional file 1: Figures S1, S2). An elongation temperature of 60 °C reduced the LOD by tenfold and/or produced darker bands for all exon 1/2 targets and for the first-round reaction of a single exon 2 target (assay 4). Using optimized elongation temperatures and under our lab conditions, the LOD of published PCR assays for *pfhrp2* and *pfhrp3* varied, with lower limits ranging from approximately 0.01 pg/ μ L to 0.1 ng/ μ L of 3D7 DNA, or approximately 0.4–4000 parasite genomes/ μ L (Table 1). Additionally, while adding a second round of amplification in the best performing nested PCR for *pfhrp2* exon 1/2 (assay 1) resulted in darker bands, the assay's LOD was unchanged.

Amplification of paralogous genes by exon 1/2 assays

Unexpectedly, paralogous amplification by assays targeting exon 1/2 but not exon 2 of both genes was observed. Bands were visualized from assays targeting exon 1/2 of *pfhrp2* (assay 1) using Dd2 (*pfhrp2*-deleted) control DNA (Additional file 1: Figure S3) and those targeting exon 1/2 of *pfhrp3* (assay 5) using HB3 (*pfhrp3*-deleted) control DNA (Additional file 1: Figure S4), respectively. Assay 1 produced unexpected bands at all three extension

Table 1 Published *pfhrp2/3* primer sequences, limits of detection, and optimized conditions

Assay #	Target	Primer sequences (5'→3')	Reaction conditions	Cycling parameters	LOD,* ng/μL (genomes/μL)	Ref
1	<i>pfhrp2</i> (exon 1/2)	Outer: For: GGTTCCTTCTCAAAAA TAAAG Rev: TCTACATGTGCTTGA GTTTCG Optional—inner: For: GTATTATCCGCTGCCGTT TTTGCC Rev: CTACACAAGTTATTATTA AATGCGGAA	200 nM each primer HotStarTaq MM (Qiagen, Venlo, Netherlands) 3 μL template DNA or 100 × diluted first-round product 25 μL reaction vol	95 °C × 15 min; 40 cycles of 94 °C × 1 min, 50 °C × 1 min for outer primers or 55 °C × 1 min for inner primers, 60 °C × 1 min; 60 °C × 10 min	10 ⁻⁵ (~0.4)	[3]
2	<i>pfhrp2</i> (exons 1/2)	For: TATCCGCTGCCGTTT TTGCC Rev: AGCATGATGGGCATC ATCCTA	400 nM each primer HotStar- Taq MM 3 μL template DNA 25 μL reaction vol	95 °C × 15 min; 40 cycles of 94 °C × 1 min, 57 °C × 1 min, 60 °C × 1 min; 60 °C × 10 min	10 ⁻¹ (~4000)	[9]
3	<i>pfhrp2</i> (exon 2)	For: ATTCCGCATTTAATAATA ACTTGTGTAGC Rev: ATGGCGTAGGCAATG TGTGG	400 nM each primer HotStarTaq MM 3 μL template DNA 25 μL reaction vol	95 °C × 15 min; 40 cycles of 94 °C × 1 min, 59 °C × 1 min, 72 °C × 1 min; 72 °C × 10 min	10 ⁻⁴ (~4)	[5]
4	<i>pfhrp2</i> (exon 2)	Outer: For: CAAAAGGACTTAATTTAA ATAAGAG Rev: AATAAATTAATGGCGTA GGCA Optional—inner (hemi- nested): For: ATTATTACACGAAACTCA AGCAC Rev: AATAAATTAATGGCGTA GGCA	400 nM each primer HotStarTaq MM 3 μL template DNA or 100× diluted first-round product 25 μL reaction vol	95 °C × 15 min; 40 cycles of 94 °C × 1 min, 57 °C × 1 min for outer primers or 62 °C for inner primers, 60 °C × 1 min; 60 °C × 10 min	10 ⁻⁴ (~4)	[23]
5	<i>pfhrp3</i> (exons 1/2)	For: TATCCGCTGCCGTTTTTG CTTCC Rev: TGCATGATGGGCATC ACCTG	400 nM primers HotStarTaq MM 3 μL template DNA 25 μL reaction vol	95 °C × 15 min; 40 cycles of 94 °C × 1 min, 60 °C × 1 min, 60 °C × 1 min; 60 °C × 10 min	10 ⁻⁴ (~4)	[9]
6	<i>pfhrp3</i> (exon 2)	Outer: For: AATGC AAAAGGACTT AATTC Rev: TGGTGTAAGTGATGC GTAGT Optional—inner (hemi- nested): For: AAATAAGAGATTATTACA CGAAAG Rev: TGGTGTAAGTGATGC GTAGT	400 nM primers HotStarTaq MM 3 μL template DNA or 100 × diluted first-round product 25 μL reaction vol	95 °C × 15 min; 40 cycles of 94 °C × 1 min, 55 °C × 1 min, 60 °C × 1 min; 60 °C × 10 min	10 ⁻⁴ (~4)	[23]
7	<i>pfldh</i> (initial qPCR)	For: ACGATTGGCTGGAG CAGAT Rev: TCTCTATTCCATTCTTTG TCACTCTTTC Probe: FAM-GTAATAGTAACA GCTGGATTACCAAGGCC CCA-TAMRA	200 nM primers 100 nM probe Probe Master qPCR Mix (Roche Diagnostics, Indian- apolis, IN) 2 μL template DNA 12 μL reaction vol	50 °C × 2 min; 95 °C × 10 min; 40 cycles of 95 °C × 15 s, 60 °C × 1 min	10 ⁻⁴ (~4)	[31]
8	<i>Pf-tubulin</i> (confirmatory)	For: AATAAATCATAATGATGT GCGCAAGTGATCC Rev: AATAAATCATAATCCTTT GTGGACATTCTCCTC	300 nM primers FastStart Universal SYBR Green MM (Roche Diagnostics) 3 μL template DNA 25 μL reaction volume	50 °C × 2 min; 95 °C × 10 min; 40 cycles of 95 °C × 15 s, 60 °C × 1 min; Dissociation analysis	10 ⁻³ (~40)	[32, 33]

LOD lower limit of detection; MM master mix; Vol volume; Bp base pair

* Typical LOD under the conditions of this laboratory. Assay performance varied between runs, but consistently achieved LODs within one log₁₀ of the listed LOD

temperatures tested, while assay 5 produced bands when tested using 60 and 65 °C extension temperatures but not at 72 °C. All tested exon 2 assays for both genes (assay 3, 4, and 6) produced negative results using Dd2 or HB3 DNA, as expected. *Pfhrp2* exon 1/2 assay 2 was not tested due to its poor performance during initial LOD testing, presumably a result of a single base insertion near the 3' end of the reverse primer compared to the 3D7 reference sequence. The sequence homology of *pfhrp2* and *pfhrp3* exon 1/2 primer binding sites (Fig. 1) suggested that amplification of paralogous genes had occurred—i.e., the amplicons generated by the *pfhrp2* exon 1/2 assay using Dd2 strain control DNA represented *pfhrp3* amplification and that those generated by the *pfhrp3* exon 1/2 assay using HB3 strain control DNA represented *pfhrp2* amplification.

Sequencing results confirmed amplification of *pfhrp3* by the *pfhrp2* exon 1/2 assay, and vice versa. With Dd2

strain (*pfhrp2*-deleted) template, the *pfhrp2* exon 1/2 assay (assay 1) unexpectedly produced a single band with a fragment length of approximately 300 bp. The amplicon's sequence aligned to the *pfhrp3* gene with 92% sequence homology (Additional file 1: Figure S5). With HB3 strain (*pfhrp3*-deleted) template, the *pfhrp3* exon 1/2 assay (assay 5) unexpectedly produced two clear bands with fragment lengths of approximately 300 and 800 bp and a faint band at approximately 400 bp. Sequences generated using DNA extracted from each band aligned to the *pfhrp2* gene, with 98, 99, and 97% sequence homology for the 300, 400, and 800 bp fragments, respectively (Additional file 1: Figure S6). However, when applied to 3D7 strain (*pfhrp2/3*-positive) control template, both exon 1/2 assays produced the expected result: a single band with a sequence that aligned to *pfhrp2* with 96% homology or *pfhrp3* with 99% sequence homology (for assays 1 and 5, respectively).




Fig. 1 *Pfhrp2* and *pfhrp3* gene sequence homology. Alignment of reference sequences from the consensus 3D7 (v3.0) genome, 5'→3', with expected binding sites for *pfhrp2* assays (white boxes) and *pfhrp3* assays (gray boxes). The reverse primer sequence for Assay 2 includes a single-base insertion (cytosine) at the location indicated by an asterisk (*). Identical bases are indicated by a period (.), missing bases by a dash (-), substitutions by the discordant base, and non-coding regions by lower case font

Streamlined testing algorithm

These findings were used to develop a streamlined, PCR-based testing pipeline for *pfhrp2/3*-negative *P. falciparum* (Fig. 2) that incorporates optimized elongation temperatures, LOD testing results, and assay specificity.

Discussion

By lowering the LOD and employing assays that distinguish *pfhrp2* from *pfhrp3*, this testing algorithm provides an improved approach to PCR-based detection of *pfhrp2/3*-negative *P. falciparum*. Importantly, PCR-based approaches for identification of *pfhrp2*- and *pfhrp3*-negative parasites must be coupled with verification of *P. falciparum* parasitaemia and confirmation that parasite DNA is present at concentrations above the LODs of the *pfhrp2* and *pfhrp3* assays. These goals were achieved by employing assays targeting two *P. falciparum*-specific, single-copy genes, lactate dehydrogenase (*pfldh*) and *P.*

falciparum beta tubulin (PfBtubulin), as the initial and final steps of the testing pipeline.

Lowering the elongation temperature improved the LOD of all published assays with exon 1/2 targets on either gene. This finding likely represents improved amplicon extension across the AT-rich intron between the exons as suggested by previous reports [22]. Unexpected amplification of paralogous gene targets by the *pfhrp2* and *pfhrp3* exon 1/2 assays was observed, presumably due to sequence homology at the primer binding sites. In regions where co-existing *pfhrp2* and *pfhrp3* deletions are common, the impact of non-specific amplification is expected to be reduced [25–27]. Additionally, the absence of paralogous amplification of 3D7 control DNA suggests that the availability of abundant, completely homologous primer binding sites early in PCR cycling reduces the likelihood of exponential amplification after mispriming. To reduce the risk of unintentional amplification of paralogous




Fig. 2 Proposed testing pipeline. Single-step assays are favoured to reduce the risk of contamination, and assays should be performed in duplicate. In addition to positive *P. falciparum* (e.g. 3D7 strain) DNA and no template controls, either **pfhrp2*-negative (e.g. Dd2 strain) or ***pfhrp3*-negative (e.g. HB3 strain) *P. falciparum* DNA controls should be used

genes, this testing algorithm uses assays targeting exon 2 of both genes. This approach also permits analysis of the repetitive sequences that encode epitopes recognized by anti-PfHRP2 antibodies [28].

A broad range of LOD results was observed for published *pfhrp2* and *pfhrp3* assays, spanning over 4 orders of magnitude under the laboratory conditions employed during this study. These differences were addressed in the resulting testing pipeline (Fig. 2) by defining an initial threshold DNA concentration ten-fold higher than the LOD of the downstream *pfhrp2* and *pfhrp3* assays. In addition, a stringent, final single-copy-gene PCR that meets the same LOD requirement was included, providing confirmation that sample degradation has not occurred during the testing process. The typical workflow employed in this laboratory includes assays 3 and 6 for *pfhrp2* and *pfhrp3* testing, respectively, performed in duplicate. For discordant results (i.e. one of two replicates positive), samples are called positive if there is a clear band of appropriate fragment length. If not, the assay is repeated, and the final call is based on the third result. Because the first-round of the nested assays achieved LODs below the initial and final confirmatory, falciparum-specific assays, their use as single-step assays is favoured in this laboratory to reduce the risk of contamination and improve work flow.

In settings where real-time PCR is not feasible, the proposed initial lactate dehydrogenase (*pfldh*) quantitative PCR assay and the final confirmatory *P. falciparum* beta tubulin (*PfBtubulin*) assays could be replaced with traditional PCR assays with LODs above the LOD of the *pfhrp2* and *pfhrp3* assays. Because assay performance can vary from laboratory to laboratory and with different reagents or equipment, it is essential to confirm the LOD of each assay using the reagents and laboratory infrastructure at hand.

In addition to PCR-based testing, current guidelines recommend independent confirmation of *P. falciparum* parasitaemia using microscopy or a non-PfHRP2-based RDT, such as an RDT that detects *P. falciparum* lactate dehydrogenase (pf-pLDH), before making deletion calls [19, 21]. Quantification of circulating PfHRP2 antigen is also a valuable tool that can be particularly useful for assessing PfHRP2-RDT-negative but *pfhrp2/3*-PCR-positive isolates with impaired protein expression [29]. Additionally, novel assays under development such as those targeting regionally specific deletion breakpoints or employing droplet digital PCR, have potential to improve throughput [30].

Conclusions

Surveillance of *pfhrp2*- and *pfhrp3*-negative *P. falciparum* requires careful laboratory workflows. PCR-based testing methods, coupled with microscopy and/or antigen testing, serve as useful tools to support policy development. Standardized approaches to the detection of *pfhrp2*- and *pfhrp3*-negative *P. falciparum* should inform efforts to define the impact of these parasites [20, 21].

Additional files

Additional file 1: Figure S1. *Pfhrp2* assay performance using serially diluted *P. falciparum* 3D7 strain DNA. Elongation temperatures were varied as listed below. All other reaction conditions are specified in Table 1.

Figure S2. *Pfhrp3* assay performance using serially diluted *P. falciparum* 3D7 strain DNA. Elongation temperatures were varied as listed below. All other reaction conditions are specified in Table 1. **Figure S3.** Representative agarose gel electrophoresis depicting unexpected spurious bands from Dd2 strain (*pfhrp2*-deleted) control DNA. PCR targeting *pfhrp2* exon 1/2 (assay 1 outer) yielded a spurious ~300 bp band from serial dilutions of *pfhrp2*-deleted Dd2 strain control DNA at all three elongation temperatures. **Figure S4.** Representative agarose gel electrophoresis depicting unexpected spurious bands from HB3 strain (*pfhrp3*-deleted) control DNA. PCR targeting *pfhrp3* exon 1/2 (assay 5) yielded spurious bands at ~300, 400, and 800 bp from serial dilutions of Dd2 strain control DNA using optimized elongation temperatures (Table 1). **Figure S5.**

Pfhrp3 assay performance using serially diluted *P. falciparum* 3D7 strain DNA. Elongation temperatures were varied as listed below. All other reaction conditions are specified in Table 1. The sequence of *pfhrp2* exon 1/2 (assay 1) PCR product aligns to *pfhrp3*, due to spurious PCR amplification of the Dd2 *pfhrp3* gene. PCR was performed using 0.01 ng/μL of 3D7 (*pfhrp2*-positive) and Dd2 (*pfhrp2*-negative) control DNA, respectively (see Additional file 1: Figure S3), followed by Sanger sequencing of amplicons. Reference sequences from the consensus 3D7 (v3.0) genome for *pfhrp2* and *pfhrp3* are displayed on the top two rows (REF), from 5'→3', with capital letters for coding regions and genetic coordinates in reference to the *pfhrp2* gene. Identical bases are indicated by a period (.), missing bases by a dash (-), substitutions by the discordant base. PCR product sequence contigs are highlighted as follows: 3D7 control DNA (light gray); Dd2 control DNA (dark gray). **Figure S6.** The sequences of *pfhrp3* exon 1/2 (assay 5) PCR product align to *pfhrp2*, due to spurious PCR amplification of the HB3 *pfhrp2* gene. PCR was performed using 0.01 ng/μL of 3D7 (*pfhrp3*-positive) and HB3 (*pfhrp3*-negative) control DNA, respectively (see Additional file 1: Figure S4), followed by Sanger sequencing of amplicons. Reference sequences from the consensus 3D7 (v3.0) genome for *pfhrp2* and *pfhrp3* are displayed on the top two rows (REF), from 5'→3', with capital letters for coding regions and genetic coordinates in reference to the *pfhrp3* gene. Identical bases are indicated by a period (.), missing bases by a dash (-), substitutions by the discordant base. PCR product sequence contigs are highlighted as follows: 3D7 control DNA (medium gray) and HB3 control DNA 300 bp fragment (light gray), 400 bp fragment (medium gray), and 800 bp fragment (dark gray)

Additional file 2. FASTA file containing Sanger sequences of amplicons produced by *pfhrp2/3* exon 1/2 assays when applied to Dd2 (*pfhrp2*-negative), HB3 (*pfhrp3*-negative), and 3D7 (*pfhrp2/3*-positive) DNA

Authors' contributions

JBP designed the study, OA and JBP performed the laboratory analyses, analysed and interpreted the data and drafted the manuscript. SRM and JJJ made substantial contributions to its conception and design. All authors read and approved the final manuscript.

Author details

¹ Division of Infectious Diseases, Department of Medicine, University of North Carolina, 130 Mason Farm Rd, Chapel Hill, NC 27599, USA. ² Department

of Epidemiology, Gillings School of Global Public Health, University of North Carolina, 135 Dauer Dr, Chapel Hill, NC 27599, USA. ³ Curriculum in Genetics and Microbiology, University of North Carolina, 321 South Columbia Street, Chapel Hill, NC 27599, USA.

Acknowledgements

The authors wish to thank Andreea Waltmann for helpful advice and Kyaw Thwai for assistance with the assays. The following reagents were obtained through BEI Resources, NIAID, NIH: Genomic DNA from *Plasmodium falciparum*, Strain HB3, MRA-155G, contributed by Thomas E. Wellem; Genomic DNA from *Plasmodium falciparum*, Strain Dd2, MRA-150G, contributed by David Walliker.

Competing interests

The authors declare that they have no competing interests.

Availability of data and materials

Representative gel images can be found in the Additional file. Upload of Sanger sequencing contigs to GenBank was not permitted due to their short < 200 bp fragment length. Sequence alignments are displayed in Additional file 1: Figures S5 and S6, and raw sequencing reads are provided in a combined FASTA file (Additional file 2). Otherwise, data sharing is not applicable to this article, as no datasets were generated or analysed during the current study.

Consent for publication

Not applicable.

Ethics approval and consent to participate

Not applicable.

Funding

This work was supported by the National Institute of Allergy and Infectious Diseases [5R01AI107949 to SRM], the American Society of Tropical Medicine and Hygiene-Burroughs Wellcome Fund to JBP, and the Thrasher Research Fund to JBP. The funders had no role in the study design, data collection and interpretation, writing, or the decision to submit the work for publication.

Publisher's Note

Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

Received: 25 October 2017 Accepted: 21 March 2018

Published online: 02 April 2018

References

- WHO. World malaria report 2016. Geneva: World Health Organization; 2016.
- Abdallah JF, Okoth SA, Fontecha GA, Torres RE, Banegas EI, Matute ML, et al. Prevalence of *pfhrp2* and *pfhrp3* gene deletions in Puerto Lempira, Honduras. *Malar J*. 2015;14:19.
- Akinyi S, Hayden T, Gamboa D, Torres K, Bendezu J, Abdallah JF, et al. Multiple genetic origins of histidine-rich protein 2 gene deletion in *Plasmodium falciparum* parasites from Peru. *Sci Rep*. 2013;3:2797.
- Baldeviano GC, Okoth SA, Arrospe N, Gonzalez RV, Sanchez JF, Macedo S, et al. Molecular epidemiology of *Plasmodium falciparum* malaria outbreak, Tumbes, Peru, 2010–2012. *Emerg Infect Dis*. 2015;21:797–803.
- Koita OA, Doumbo OK, Ouattara A, Tall LK, Konare A, Diakite M, et al. False-negative rapid diagnostic tests for malaria and deletion of the histidine-rich repeat region of the *hrp2* gene. *Am J Trop Med Hyg*. 2012;86:194–8.
- Murillo Solano C, Akinyi Okoth S, Abdallah JF, Pava Z, Dorado E, Incardona S, et al. Deletion of *Plasmodium falciparum* histidine-rich protein 2 (*pfhrp2*) and histidine-rich protein 3 (*pfhrp3*) genes in Colombian parasites. *PLoS ONE*. 2015;10:e0131576.
- Wurtz N, Fall B, Bui K, Pascual A, Fall M, Camara C, et al. *Pfhrp2* and *pfhrp3* polymorphisms in *Plasmodium falciparum* isolates from Dakar, Senegal: impact on rapid malaria diagnostic tests. *Malar J*. 2013;12:34.
- Kumar N, Pande V, Bhatt RM, Shah NK, Mishra N, Srivastava B, et al. Genetic deletion of HRP2 and HRP3 in Indian *Plasmodium falciparum* population and false negative malaria rapid diagnostic test. *Acta Trop*. 2013;125:119–21.
- Gamboa D, Ho MF, Bendezu J, Torres K, Chiodini PL, Barnwell JW, et al. A large proportion of *P. falciparum* isolates in the Amazon region of Peru lack *pfhrp2* and *pfhrp3*: implications for malaria rapid diagnostic tests. *PLoS ONE*. 2010;5:e8091.
- Li P, Xing H, Zhao Z, Yang Z, Cao Y, Li W, et al. Genetic diversity of *Plasmodium falciparum* histidine-rich protein 2 in the China-Myanmar border area. *Acta Trop*. 2015;152:26–31.
- Amoah LE, Abankwa J, Oppong A. *Plasmodium falciparum* histidine rich protein-2 diversity and the implications for PfHRP 2: based malaria rapid diagnostic tests in Ghana. *Malar J*. 2016;15:101.
- Bharti PK, Chandel HS, Ahmad A, Krishna S, Udhayakumar V, Singh N. Prevalence of *pfhrp2* and/or *pfhrp3* gene deletion in *Plasmodium falciparum* population in eight highly endemic states in India. *PLoS ONE*. 2016;11:e0157949.
- Berhane A, Russom M, Bahta I, Hagos F, Ghirmai M, Uqubay S. Rapid diagnostic tests failing to detect *Plasmodium falciparum* infections in Eritrea: an investigation of reported false negative RDT results. *Malar J*. 2017;16:105.
- Parr JB, Verity R, Doctor SM, Janko M, Carey-Ewend K, Turman BJ, et al. *Pfhrp2*-Deleted *Plasmodium falciparum* parasites in the Democratic Republic of the Congo: a national cross-sectional survey. *J Infect Dis*. 2017;216:36–44.
- Rachid Viana GM, Akinyi Okoth S, Silva-Flannery L, Lima Barbosa DR, Macedo de Oliveira A, Goldman IF, et al. Histidine-rich protein 2 (*pfhrp2*) and *pfhrp3* gene deletions in *Plasmodium falciparum* isolates from select sites in Brazil and Bolivia. *PLoS ONE*. 2017;12:e0171150.
- Beshir KB, Sepulveda N, Bharmal J, Robinson A, Mwanguzi J, Busula AO, et al. *Plasmodium falciparum* parasites with histidine-rich protein 2 (*pfhrp2*) and *pfhrp3* gene deletions in two endemic regions of Kenya. *Sci Rep*. 2017;7:14718.
- Kozycki CT, Umulisa N, Rulisa S, Mwikarago EI, Musabyimana JP, Habimana JP, et al. False-negative malaria rapid diagnostic tests in Rwanda: impact of *Plasmodium falciparum* isolates lacking *hrp2* and declining malaria transmission. *Malar J*. 2017;16:123.
- Gupta H, Matambisso G, Galatas B, Cistero P, Nhamussua L, Simone W, et al. Molecular surveillance of *pfhrp2* and *pfhrp3* deletions in *Plasmodium falciparum* isolates from Mozambique. *Malar J*. 2017;16:416.
- WHO. *P. falciparum* *hrp2/3* gene deletions: conclusions and recommendations of a technical consultation. Geneva: World Health Organization; 2016.
- WHO. False-negative RDT results and implications of new *P. falciparum* histidine-rich protein 2/3 gene deletions. Geneva: World Health Organization; 2016.
- Cheng Q, Gatton ML, Barnwell J, Chiodini P, McCarthy J, Bell D, et al. *Plasmodium falciparum* parasites lacking histidine-rich protein 2 and 3: a review and recommendations for accurate reporting. *Malar J*. 2014;13:283.
- Su XZ, Wu Y, Sifri CD, Wellem TE. Reduced extension temperatures required for PCR amplification of extremely A + T-rich DNA. *Nucleic Acids Res*. 1996;24:1574–5.
- Baker J, McCarthy J, Gatton M, Kyle DE, Belizario V, Luchavez J, et al. Genetic diversity of *Plasmodium falciparum* histidine-rich protein 2 (PfHRP2) and its effect on the performance of PfHRP2-based rapid diagnostic tests. *J Infect Dis*. 2005;192:870–7.
- Li W, Cowley A, Uludag M, Gur T, McWilliam H, Squizzato S, et al. The EMBL-EBI bioinformatics web and programmatic tools framework. *Nucleic Acids Res*. 2015;43:W580–4.
- Berhane A, Mihreteab S, Mohammed S, Embaye G, Hagos F, Zehaie A, et al. PfHRP2 detecting malaria RDTs: alarming false-negative results in Eritrea (Poster #879). Annual meeting of the American Society of Tropical Medicine and Hygiene. 2016; Atlanta, Georgia, USA.
- Gamboa D, Ho MF, Bendezu J, Torres K, Chiodini PL, Barnwell JW, et al. A large proportion of *P. falciparum* isolates in the Amazon region of Peru lack *pfhrp2* and *pfhrp3*: implications for malaria rapid diagnostic tests. *PLoS ONE*. 2010;5:e008091.

27. Menegon M, L'Episcopia M, Nurahmed AM, Talha AA, Nour BYM, Severini C. Identification of *Plasmodium falciparum* isolates lacking histidine-rich protein 2 and 3 in Eritrea. *Infect Genet Evol*. 2017;55:131–4.
28. Lee N, Gatton ML, Pelecanos A, Bubb M, Gonzalez I, Bell D, et al. Identification of optimal epitopes for *Plasmodium falciparum* rapid diagnostic tests that target histidine-rich proteins 2 and 3. *J Clin Microbiol*. 2012;50:1397–405.
29. Rogier E, Plucinski M, Lucchi N, Mace K, Chang M, Lemoine JF, et al. Bead-based immunoassay allows sub-picogram detection of histidine-rich protein 2 from *Plasmodium falciparum* and estimates reliability of malaria rapid diagnostic tests. *PLoS ONE*. 2017;12:e0172139.
30. Koepfli C, Badu K, Lo E, Hemming-Schroeder E, Yan G. High-throughput screening of *P. falciparum* *hrp2* deletion for monitoring of RDT efficacy for malaria diagnosis (Poster #290). Annual meeting of the American Society of Tropical Medicine and Hygiene. 2017; Baltimore, Maryland, USA.
31. Pickard AL, Wongsrichanalai C, Purfield A, Kamwendo D, Emery K, Zalewski C, Kawamoto F, Miller RS, Meshnick SR. Resistance to antimalarials in Southeast Asia and genetic polymorphisms in *pfmdr1*. *Antimicrob Agents Chemother*. 2003;47:2418–23.
32. Price RN, Uhlemann AC, Brockman A, McGready R, Ashley E, Phaipun L, et al. Mefloquine resistance in *Plasmodium falciparum* and increased *pfmdr1* gene copy number. *Lancet*. 2004;364:438–47.
33. Afonina I, Ankoudinova I, Mills A, Lokhov S, Huynh P, Mahoney W. Primers with 5' flaps improve real-time PCR. *Biotechniques*. 2007;43:770–4.

Submit your next manuscript to BioMed Central
and we will help you at every step:

- We accept pre-submission inquiries
- Our selector tool helps you to find the most relevant journal
- We provide round the clock customer support
- Convenient online submission
- Thorough peer review
- Inclusion in PubMed and all major indexing services
- Maximum visibility for your research

Submit your manuscript at
www.biomedcentral.com/submit

